

## GROWING BACTERIA CULTURES: GENERAL

1. L Broth prep:
  - There are many different broths for growing bacteria. For most protocols performed in this laboratory we use L broth (LB) as media for growing bacteria cells in solution and/or on plates
  - We use LB tablets
  - 1 tablet / 50ml dH<sub>2</sub>O
  - In a Erlenmeyer flask mix water and tablets
  - Cover top with a layer of tin foil
  - Autoclave flask
  - Cool broth before addition of antibiotics
  
2. Agar Plates for growing bacteria
  - 7.5g Agar
  - 500ml LB ( or 500ml H<sub>2</sub>O + 10 tablets LB)
  - Mix together in an Erlenmeyer
  - Tin foil top and autoclave
  - Let cool to about 50°C then add antibiotics
  - Pour agar into Petri dishes just to cover bottom of plate
  
3. Addition of Antibiotics
  - In this lab we mostly use Ampicillin [0.1mg/ml]
  - The stock solutions of Amp [25mg/ml]
  - Example of calculation
$$300\text{ml culture: } C_1V_1 = C_2V_2$$
$$(300\text{ml})(0.1\text{mg/ml}) = (X\text{ml})(25\text{mg/ml})$$
$$X = 1.2\text{mls amp in } 300\text{ml LB}$$
  
4. Overnight Cultures (for plasmid preps)
  - Add 5ml of LB broth to a 15ml conical tube.
  - Add 20ul ampicillin [stock 25mg/ml] [final 100ug/ml]
  - Pick a colony from the plate with loop and shake loop in LB amp broth or add 200µl bacterial glycerol.
  - Repeat for all plates that DNA will be obtained from.
  - Place all the tubes in the warm room and shake at 250rpm's overnight.
  
5. Glycerols
  - Inoculate 5ml of LBroth (w/o antibiotics) with 300µl of an overnight culture
  - Incubate with shaking at 37°C room for about 1-1.5 hour (mid-log phase)
  - Add equal amount of sterile glycerol solution (5ml) invert mix.
    - 65% glycerol
    - 0.1M MgSO<sub>4</sub>
    - 0.025M Tris-Cl (pH 8)
  - Aliquot 1ml of bacteria glycerol into a cryovials
  - Store -80°C

NB: make 10 glycerols for each new plasmid generated. This stock is used for inoculating for plasmid preps