

IP PROTOCOL USING STREPAVIDIN AGAROSE

NB: Remember to keep everything RNase and DNase free, as you will be extracting RNA.

1. Make HNTM Buffer (500ml):

	MW	[final]	Amt
HEPES pH 7.5	238.3	50mM	5.95g
NaCl	58.44	150mM	4.383g
Triton X-100		1%	5ml
MgCl₂		1mM	0.5ml

*Also make 1L of HNM buffer (without Triton X-100)

2. Use the BCA (Pierce) Protein Assay to determine the concentration of the cytosolic and nuclear lysate. Calculate the volume of lysate needed, using 75 μ g of cytosolic lysate and 37.5 μ g of nuclear lysate to immunoprecipitate. Typically we use 6 tubes to take to RNA and 1 tube as a western control. You will need to bring the volume to 850 μ l using HNM buffer.
3. The Streptavidin Agarose comes in a 50/50 slurry, but not in the correct buffer. To change the buffer of the beads, determine the quantity that will be used for the samples that will be immunoprecipitated (100ul of 50/50 slurry per tube). Spin the correct volume needed down at 4°C and then aspirate off the supernatant. Resuspend the bead with a generous amount of HNM buffer, spin down, and aspirate off the supernatant again. This is repeated 3-5 times. Any sodium azide that was in the original buffer that may still be present will be washed with HNTM buffer. Wash as above with HNTM buffer 2-3 times. Then add an equal amount of HNM buffer to the beads to acquire a 50/50 slurry in HNM buffer.
4. For each tube, add 100ul of the 1:1 bead mixture and spin down at 13,200 rpm at 4°C. Make sure there is approximately 50ul of beads per tube.
5. Add the appropriate volume of HNM buffer needed for each sample of beads to reach 850ul of HNM/lysate, and then add the following to an RNase/DNase free microfuge tube in this order:

Halt Protease Inhibitor (Pierce Cat # 78410)	10 μ l
Vanadyl Ribonucleoside Complex	10 μ l
Rnase Out 40U/ μ l (Invitrogen Cat # 10777-019)	10 μ l
0.1M DTT	1 μ l
0.5M EDTA	33 μ l

6. Now add the appropriate volume of lysate based on your calculation for the final concentration.
7. Rotate for 3 hours at room temperature.
8. Spin down the bead with lysate (13,200rpm at 4C) and wash 3 times with cold 1M Urea (MW=60.06g) in HNTM buffer. Vortex the sample and Urea and let sit 1-2 minutes before spinning again.
9. Wash beads 2 times with cold HNTM buffer.
10. Take a sample for western control.