

POLYMERASE CHAIN REACTION USING TAO BEADS

DNA template: Use cDNA from Reverse Transcription reaction or plasmid DNA for amplification in PCR.

Reagents Needed:

Forward and reverse primers (10 μ M)
puReTaq Ready-To-Go PCR beads in 0.5ml tubes (Amersham Biosciences Cat # 27-9558-01)

Protocol

1. Resuspend the PCR bead in 18 μ l of RNase/ DNase free H₂O. Let bead dissolve.
2. Add 5 μ l of cDNA from RT , If using plasmid DNA, add 100ng (if amount exceeds 5 μ l, resuspend bead in less dH₂O, the object being to have a 25 μ l total volume reaction)
3. May/ may not need to include positive control (ask)
4. Negative control : Resuspend bead in 23 μ l dH₂O (no DNA)
5. Add 1 μ l (10 μ M) each of **forward** and **reverse** primers.
6. Vortex gently and spin down the contents of the tube.
7. Cycle on the PCR program on the thermal cycler in 5S-19 (PCR program): program is preset except for entering your particular melting temp (step c)...can be setup to run overnight.
 - a. 80°C for 5 minutes
 - b. 94°C for 1 minute
 - c. X°C for 1 minute (use temp 5 degrees lower than lowest melting temp of your
 - d. 72°C for 1 minute
 - e. Repeat b-d for 27 cycles
 - f. 72°C for 10 minutes
 - g. 4°C ∞

DNA Gel:

1. Pour a 1% agarose gel w/ ethidium bromide in gel (5 μ l/ 100 ml gel; 10 mg/ml EtBr stock) with a 10-lane comb (typical). Use 1x TAE buffer for reservoir.
2. Add 5 μ l of **DNA LB WITHOUT dye** to each PCR product vial. Don't heat samples!
3. Load your samples on gel. Include wells with
 - a. 8 μ l of Invitrogen 100bp ladder prepared solution
 - b. Dye only control well with Loading buffer WITH dye (20 μ l water + 5 μ l LB w/ dye)
4. Run gel at 100 V for ~90-120 min.

Take a photo of your gel using camera in Dr. Shewchuck's lab.
Determine size of PCR product using the ladder.