ISOLATION OF TOTAL RNA FROM 10 CM² PLATES

Materials: TRIzol Reagent InVitrogen cat#15596-026 (100mL) Scrapers (sterile, RNase free), Syringes and 20g needles, sterile 15ml polypropylene tubes or sterile microfuge tubes.

Cell Culture:

- (3) 10cm plates for Preadipocyte time points
- (2) 10cm plates for early differentiation time points
- (1) 10cm plates for Adipocytes

The Cell Harvest:

- 1. Aspirate media from plates.
- 2. Wash plates once with ice cold 1x PBS and aspirate off PBS
- 3. Add 1ml TRIzol reagent per plate.
- 4. Scrape cells with cell scraper and pool cell mass into corner of plate
- 5. Triturate cell suspension with 20g needle on a 5ml sterile syringe.
- 6. The lysate is then transferred to the second plate and the procedure repeated. And so forth with the third plate.

Example: 3 ml of TRIzol on the first plate, scraping the cell mass to one corner, then pipetting it up and down using a sterile syringe with a 20-gauge needle. This insures lysis of the cells and shears the DNA. Then transfer to next plate.

7. Resulting suspension is 3 plates in 3 ml, which is then transferred to a 15 ml polypropylene tube and stored at -80° until the time of RNA isolation.

RNA Isolation

- 1. Thaw homogenates and aliquot out 1ml into 1.5ml microfuge tubes
- 2. Incubate at room temp for 5minutes
- 3. Add 200ul Chloroform, vortex into single phase & incubate for 3min @RT
- 4. Centrifuge 12k x g for 15min @ 4°C
- 5. Transfer aqueous phase (top layer) to clean tube
- 6. Optional: second chloroform extraction to remove residual phenol
- 7. Add 500ul Isopropanol and vortex
- 8. Incubate 10 min at room temp (can be froze -80°C o/n) precip. often visible
- 9. Centrifuge 12k x g for 10min @ 4°C

- 10. Remove supernatant
- 11. Wash with 1ml 75% ethanol (add to tube and gently vortex)
- 12. Centrifuge 7.5k x g 5min @ 4°C
- 13. Aspirate off supernatant
- 14. Air dry pellet (pellet is white while wet and dries clear, don't over dry)
- 15. Resuspended pellet in 70ul of nuclease free dH₂O
- 16. Incubate 10 min @ 55°C
- 17. Place on ice
- 18. NanoDrop for concentration
- 19. Store -80°C